

Single-cell clarity and heterogeneity in copy number profiles in primary synovial & Ewing sarcoma with ResolveDNA™ genomic amplification

Katherine A. Kennedy¹, Jon S. Zawistowski¹, Viren R. Amin¹, Durga M. Arvapalli¹, Isai Salas-González ¹, Natalie Jäger², Hanno Glimm³, Gary L. Harton¹, Stefan Fröhling², Jason A.A. West¹, Stefan Pfister².
1. BioSkryb Genomics, Inc., Durham, NC., 2. Hopp Children's Cancer Center Heidelberg (KITZ) and German Cancer Consortium (DKTK), Heidelberg, Germany, 3. NCT Dresden and German Cancer Consortium (DKFZ), Dresden, Germany

Abstract

Soft tissue sarcomas including Ewing (ES) & synovial (SS) sarcoma represent a diverse set of mesenchymal malignancies frequently driven by translocation events. While well characterized, limited data exist surrounding broader genomic alterations & intra-tumoral heterogeneity. Unfortunately, this precludes personalized prognoses and treatment, with 50-70% of cases resulting in relapse or progression, leading to overall poor clinical outcomes¹⁻³.

In a collaboration with the Hopp Children's Cancer Center Heidelberg (KITZ), we utilized the ResolveDNA™ process to assess copy number aberration (CNA) at the single cell level in two soft-tissue sarcoma samples, a putative Ewing sarcoma sample, and a putative synovial sarcoma.

ResolveDNA™ chemistry attenuates amplicon size, redirects amplification to the primary DNA template and avoids exponential copying of amplicons⁴. This results in unprecedented coverage and uniformity, with high SNV precision and sensitivity, allelic balance, and the enhanced ability to accurately call CNA⁴.

Single cell analysis in this study demonstrated diverse phenotypes – SS cells harboring highly disrupted genomes, with stark focal, sub-chromosomal gains and losses, and ES cells maintaining a more typical profile, with limited discrete, full-chromosome (CHR) gains and losses.

Whole exome sequencing demonstrated profile differences between the two sample types.

Methods

FACS Sorting





Patient derived cells were obtained and sorted in collaboration with the Hopp Children's Cancer Center Heidelberg. Single cells were deposited into individual wells containing 3uL of BioSkryb Genomic's proprietary cell buffer.

ResolveDNA™

Plates underwent ResolveDNA™ per manufacturer's instructions. After amplification, libraries were prepared using ResolveDNA™ library amplification & quantified using high sensitivity Qubit and an Agilent tape station. DNA libraries were equimolar pooled & sequenced on an Illumina MiniSeq targeting 2M reads per sample using 2x75 paired end sequencing chemistry. Exome samples were generated using the IDTV2 xGen panel with Kapa fragmentase after ResolveDNA™ library preparation with 12plex multiplexing. A total of 500ng library was input to the fragmentation reaction. All exomes were sequenced targeting 40M reads per sample using 2x150 paired end sequencing chemistry on an Illumina NovaSeq.

Sequencing & Data analysis

Fastq files were analyzed using multiple BaseJumper modules – including CNA and SNV assessments, as well as genome assembly and quality evaluation.



Results

Single cell analysis of ES demonstrated significant clonality with minimal CNA. All cells harbor gains on CHRs 2 and 8, & losses of 16. Greater than 50% of cells demonstrate loss of CHR 19. Unlike synovial sarcoma, minimal numbers of cells showed additional alterations. Those that did exist were focal gains on CHR 10 (n=4/12 cells), and losses on CHRs 9 and 11 (n=1 cell each).

Strikingly, copy number aberrations in the synovial sarcoma cells were numerous, with focal points of large-scale amplification and losses, including at the *MYC* locus. Intriguingly, single-cell heterogeneity was notable in multiple sub-chromosomal regions including 3q, 4p and 16q each seen in at least two independent cells. The expected SS18:SSX translocation was not identified in this patient, suggesting that these cells may have originated from a secondary sarcoma as opposed to the originally diagnosed synovial sarcoma.

Summary

Sarcomas represent a group of mesenchymal tumors frequently driven by translocation events¹. They can be difficult to treat, often resulting in recurrence and metastasis¹⁻³. Typically, these tumors have been thought of as being highly clonal, with minimal disruption of the genome. The ability to identify rare instances of sub-clonal populations of tumor cells is likely to improve prognostic outcomes in these diseases.

The heterogeneity observed in this study demonstrates the power of single cell vs bulk analysis. Here we reveal the nature of genomic instability and aid in the characterization of sarcoma subtypes at unprecedented levels of resolution. Despite the anticipated clonality of ES cells, evidence of sub-clonal heterogeneity can be seen.

The ability to detect these rare events are particularly crucial to understand disease development and progression in orphan sarcomas, especially those with complex genome rearrangements and/or poor clinical outcome. Further studies, including transcriptomic evaluation are ongoing.

References & Acknowledgments

We thank the patients & families whose donations made this research possible as well as the technical staff at the KITZ hospital.

- 1) Gazendam AM, Popovic S, Munir S, Parasi N, Wilson D, Ghert M. Synovial Sarcoma: A Clinical Review. *Curr Oncol*. 2021 May 19;28(3):1909-1920. doi: 10.3390/curronc28030177. PMID: 34009748. PMCID: PMC8161765.
- 2) Connolly EA, Bhaddi VA, Wake J, Ingley KM, Lewin J, Bae S, Wong DD, Long AP, Pryor D, Thompson SR, Strach MC, Grimison PS, Mahar A, Bonar F, Maclean F, Hong A. Systemic treatments and outcomes in CC-rearranged Sarcoma: A national multi-center clinicopathological series and literature review. *Cancer Med*. 2022 Apr;11(8):1805-1816. doi: 10.1002/cam4.4580. Epub 2022 Feb 17. PMID: 35178869. PMCID: PMC9041083.
- 3) Ffrench M, Samir A, Sponzo P, Zucchi R, Giannini C, Caldar E, Pirini MG, De Paolo M. The Biology of Synovial Sarcoma: State-of-the-Art and Future Perspectives. *Curr Treat Options Oncol*. 2021 Oct 13;22(12):109. doi: 10.1007/s11864-021-00914-4. PMID: 34687866. PMCID: PMC8541077.
- 4) González-Pérez V, Natarajan S, Xia Y, Klein D, Carter R, Pang S, Shiner B, Annu K, Putnam D, Chen W, Connolly J, Pruett-Miller S, Chen X, Easton J, Gawad C. Accurate genomic variant detection in single cells with primary template-directed amplification. *Proc Natl Acad Sci U S A*. 2023 Jun 15;120(24):e2024176118. doi: 10.1073/pnas.2024176118. PMID: 36999468. PMCID: PMC1214697.

